

eLUG Changes
9/23/2024

MAJOR CHANGES NEW TESTS

PI Analyst	TEST CODE	TEST NAME
Marilyn	WSMNCS	Spinal Muscular Atrophy Carrier Screening
Marilyn	SMNDX	Spinal Muscular Atrophy, Diagnostic Assay
Marilyn	WHSTNT	Troponin T, High Sensitivity (hs-TnT)
Marilyn	WWLNCP	Allergen- Walnut Component Panel, IgE

MAJOR CHANGES DELETE

PI Analyst	TEST CODE	TEST NAME
Pratik	WNTXS	N-TELOPEPTIDE, CROSS-LINKED, SERUM Replaced by:
Marilyn	DNA2100030	Spinal Muscular Atrophy (SMA) Carrier Screen Replaced by: WSMNCS (Spinal Muscular Atrophy Carrier Screening)
Marilyn	WTROPT	Troponin T Replaced by: WHSTNT (Troponin T, High Sensitivity (hs-TnT))
Pratik	SREQ	IA-2 Antibody Replaced by: WIA2
		Replaced by:

MINOR CHANGES (updates)

PI Analyst	TEST CODE	TEST NAME/DESCRIPTION
Marilyn	WC2729	CA 27.29 Description of change: Updated preferred container type from GO to GOS. Alternat container type updated from RA to RAS.
Marilyn	WDPHF	Phenytoin, Free Description of change: Updated preferred container type from RA to RAS. All alternate container types removed.
Marilyn	WLPA	Lipoprotein (a) Description of change: Updated preferred container type from GO to RAS. Alternate container type update: LH, EAS, EHS, or BA.

eLUG Changes
9/23/2024

Marilyn/Jackie	WPAN	Pnrma Prntl Tst w/No Mrcrdltn Pnl Description of change: Updated # 2 under special instructions: Completed Panorama questionnaire, located HERE or available to print at the bottom of the electronic Epic order.
Marilyn/Jackie	WPAN22	Pnrma Prntl Tst w 22Q11 Mrcrdltn Pnl Description of change: Updated # 2 under special instructions: Completed Panorama questionnaire, located HERE or available to print at the bottom of the electronic Epic order.
Marilyn/Jackie	WPANFP	Pnrma Prntl Tst w Extnded Mrcrdltn Pnl Description of change: Updated # 2 under special instructions: Completed Panorama questionnaire, located HERE or available to print at the bottom of the electronic Epic order.
Marilyn	WHEPCR	Hep E Virus Quant PCR Description of change: Add special instructions: Approved for gastroenterology, hepatology and liver transplant. All others require SME approval.
Marilyn	SPE2	Serum Protein Electrophoresis Description of change: Updated collect volume to 4 ml and updated minimum volume to 2 ml.
Marilyn	TROPI	Troponin I Description of change: Updated storage instructions: Centrifuge and refrigerate.
Pratik	WIA2	IA-2 Antibody Description of change: Update test name
Pratik	MPSS	Monoclonal Protein Serum Description of change: Update Specimen Container
		Description of change:
		Description of change:
		Description of change:
		Description of change: